



16900-08

Digital Titrator

Model 16900

TABLE OF CONTENTS

SPECIFICATIONS	7
OPERATION	9
GENERAL DESCRIPTION	11
1.1 Introduction	11
1.1.1 Following a Procedure for the First Time	12
1.2 Step-By-Step	13
1.3 Helpful Hints.....	19
1.3.1 To Reuse a Partially Emptied Cartridge	19
1.3.2 To Calculate Titrant Volume Used.....	19
1.3.3 To Fill Your Own Titration Cartridges	19
1.3.4 Verifying Technique.....	20
1.4 Adapting a Buret Titration to the Digital Titrator	23
1.5 Using PermaChem® Powder Pillows	24
1.6 Safety.....	25
TITRATION PROCEDURES	29
ACID-BASE	
Acid Determination	31
Base Determination	33
ACIDITY	
Methyl Orange Method	39
Phenolphthalein (Total) Method	41
ALKALINITY	
Phenolphthalein and Total Method	45
AMMONIA, HIGH RANGE (Ammonium Hydroxide)	
Ammonia Titration Procedure	55
CARBON DIOXIDE	
Using Sodium Hydroxide	61
CHELANT, FREE	
Using Magnesium Chloride.....	65
CHELANT, TOTAL	
Using Bismuth Nitrate	69

TABLE OF CONTENTS, continued

CHLORIDE

Mercuric Nitrate Method	73
Silver Nitrate Method	74

CHLORINE, FREE AND TOTAL

DPD-FEAS Method.....	81
----------------------	----

CHLORINE, TOTAL

Iodometric Method (1 to 400 mg/L as Cl ₂ Using Sodium Thiosulfate).....	87
Iodometric Method (20 to 70,000 mg/L as Cl ₂ Using Sodium Thiosulfate).....	90

CHLORINE, FREE

Amperometric Forward Titration	95
--------------------------------------	----

CHLORINE, TOTAL

Amperometric Back Titration	103
-----------------------------------	-----

CHLORINE, TOTAL

Amperometric Forward Titration	115
--------------------------------------	-----

CHROMATE

Using Sodium Thiosulfate.....	123
-------------------------------	-----

HARDNESS DECISION TREE.....	129
-----------------------------	-----

HARDNESS, CALCIUM

Using EDTA	131
------------------	-----

HARDNESS, TOTAL

Using EDTA	139
------------------	-----

HARDNESS, TOTAL, SEQUENTIAL

Sequential Titration Procedure (Limited Sample)	147
---	-----

HYPOCHLORITE (Bleach)

Iodometric Method	155
-------------------------	-----

IRON

Using the TitraVer® Titration Cartridge	159
---	-----

NITRITE

Using Ceric Standard Solution.....	163
------------------------------------	-----

TABLE OF CONTENTS, continued

OXYGEN, DISSOLVED

Azide Modification of Winkler Method.....	167
Using a 300-mL BOD Bottle	167
Using a 60-mL BOD Bottle	170

SALINITY

Using Mercuric Nitrate	175
------------------------------	-----

SULFITE

Using Iodate-Iodide.....	177
--------------------------	-----

TURBIDITY STANDARDS

Preparing Turbidity-Free Water.....	181
-------------------------------------	-----

VOLATILE ACIDS

Using Sodium Hydroxide	187
------------------------------	-----

APPENDIX A

ACCURACY CHECK AND STANDARD ADDITIONS	191
---	-----

GENERAL INFORMATION	205
----------------------------------	------------

REPLACEMENT PARTS AND ACCESSORIES.....	207
--	-----

HOW TO ORDER	209
--------------------	-----

REPAIR SERVICE	211
----------------------	-----

WARRANTY	212
----------------	-----

SPECIFICATIONS

Digital Titrator

Delivery: 800 digits/mL or 0.00125 mL/digit

Accuracy*: $\pm 1\%$ for readings over 100 digits. (Uncertainty of readings is 1 digit. Most samples require more than 100 digits.)

Weight: 132 g (4.7 oz.)

Cartridges for the Digital Titrator

Volume: 13 mL

Number of tests: Most reagents are formulated to provide 100 typical titrations; the number may vary depending on sample concentration.

Weight (full): 56.75 g (2 oz.)

* Overall method accuracy includes, in addition to the Digital Titrator, other sources of error controlled by the analyst. The other sources of error include: sampling, sample volume, dilution (if required), end point detection, reagent quality, and interferences.



OPERATION

DANGER

Handling chemical samples, standards, and reagents can be dangerous. Review the necessary Material Safety Data Sheets and become familiar with all safety procedures before handling any chemicals.

DANGER

La manipulation des échantillons chimiques, étalons et réactifs peut être dangereuse. Lire les Fiches de Données de Sécurité des Produits (FDSP) et se familiariser avec toutes les procédures de sécurité avant de manipuler tous les produits chimiques.

PELIGRO

La manipulación de muestras químicas, estándares y reactivos puede ser peligrosa. Revise las fichas de seguridad de materiales y familiarícese con los procedimientos de seguridad antes de manipular productos químicos.

GEFAHR

Das Arbeiten mit chemischen Proben, Standards und Reagenzien ist mit Gefahren verbunden. Es wird dem Benutzer dieser Produkte empfohlen, sich vor der Arbeit mit sicheren Verfahrensweisen und dem richtigen Gebrauch der Chemikalien vertraut zu machen und alle entsprechenden Material Sicherheitsdatenblätter aufmerksam zu lesen.

PERIGO

A manipulação de amostras, padrões e reagentes químicos pode ser perigosa. Reveja a folha dos dados de segurança do material e familiarize-se com todos os procedimentos de segurança antes de manipular quaisquer produtos químicos.

GENERAL DESCRIPTION

1.1 Introduction


Hach's Digital Titrator is a new concept in titrimetric analysis. It is a precision dispensing device fitted with compact cartridges that contain concentrated titrants. Accurate titrations are made without the bulk and fragility of conventional burets.

A main drive screw in the Digital Titrator controls a plunger which forces the concentrated titrant from a titration cartridge in a carefully regulated flow. The titrator body is constructed of precision-molded, heavy-duty, chemical- and impact-resistant acetal plastic. Accuracy is rated at $\pm 1\%$ or better for a titration of more than 100 digits. For titrations less than 100, accuracy is ± 1 digit.

Titration solutions (titrants) are packaged in disposable polypropylene or Kynar[®] containers with Teflon-covered neoprene seals and polyethylene resealable closures to cover the cartridge tips. Each cartridge contains approximately 13 mL of titrating solution, enough for 50–100 average titrations. Titrant solutions are typically controlled to $\pm 0.5\%$ concentration with normality and tolerances listed on the label. Titrant concentrations are designed for titrations of 10 to 40 turns (100 to 400 digits) of the delivery knob. For the most commonly used concentration ranges, the digits appearing in the counter window correspond to the sample concentration.

GENERAL DESCRIPTION, continued

Figure 1 Hach Digital Titrator



Both portable and fixed-position titrations are possible with the Digital Titrator. The instrument has a grip for hand-held operation or it can be clamped to a TitraStir® Stir Plate or laboratory stand for stationary setups. See *Figure 1*.

Each Digital Titrator comes with five delivery tubes and a methods manual, which covers the most commonly tested parameters and the corresponding titrant cartridges. Right-angle (ninety-degree) delivery tubes for stationary setups are available as an optional accessory.

1.1.1 Following a Procedure for the First Time

Each method is divided into five sections: Procedure, Accuracy Check, Interferences, Summary of Method, and Reagents and Apparatus. For more information about how to select a procedure or for answers to chemical questions, see Hach's *Water Analysis Handbook* (literature 8376). For more information about chlorine measurement, also see the technical booklet titled, *Current Technology of Chlorine Analysis for Water and Wastewater* (literature 7019).

The **Procedure** details how to perform the method step-by-step. To select the appropriate sample volume and titration cartridge based on expected sample concentration, use the

GENERAL DESCRIPTION, continued

tables provided in each procedure. If the expected sample concentration is not known, start with one of the smaller sample volumes and determine its approximate concentration. Retest with the appropriate sample size.

The ranges in the table overlap to offer more flexibility. In most procedures, the number of digits used for each concentration range will be 100 to 400 digits.

To determine the actual concentration of the sample, use the correct digit multiplier for the sample volume and titration cartridge used.

Throughout the procedure, the notes will provide additional information.

The **Accuracy Check** provides a way to verify the results and determine if interferences are present. It also provides a method for checking the performance of reagents, the Digital Titrator and the operator's technique. Further information is provided in *Appendix A, Accuracy Check and Standard Additions*.

The **Interferences** section identifies common interferences causing inaccurate results and describes how to eliminate their effects. The interference levels are based on the sample volume that has 1.0 as the digit multiplier. Higher interference levels may be tolerated if a smaller sample is used.

The **Summary of Method** section discusses the chemical reaction taking place and information that applies to the entire procedure.

The **Reagents and Apparatus** list concludes the procedure. All the items required to perform the test are listed first and are available from Hach. The items listed in the notes or interferences sections are included in the optional listings.

1.2 Step-By-Step

1. Select a sample volume and titration cartridge corresponding to the expected sample concentration from the table given in each procedure.

If the expected sample concentration is not known, start with one of the smaller sample volumes and determine

GENERAL DESCRIPTION, continued

its approximate concentration. Retest with the appropriate sample size.

2. Slide the cartridge into the titrator receptacle and lock in position with a slight turn. See *Figure 2*.


Figure 2 Sliding the Cartridge into Place



3. Remove the polyethylene cap and insert a clean delivery tube into the end of the cartridge until it is tight. See *Figure 3*. Use a straight tube with a hook at the end for hand-held titrations; use a 90° tube with a hook at the end for stationary setups.

Do not insert tube past cartridge extension; see illustration below. In some instances, it might be necessary to remove a small burr on the leading edge of the tube before insertion.

Figure 3 Inserting the Delivery Tube



4. For stationary titrations, use a TitraStir Stir Plate or a clamp holder and clamp to attach the titrator to a laboratory stand. See *Figure 4* and *Figure 5*.

The TitraStir Stir Plate holds the Digital Titrator during the titration and also stirs the sample at a constant speed, leaving the analyst free to detect the end point.

GENERAL DESCRIPTION, continued

When a TitraStir Stir Plate is used, substitute or add the following Optional Apparatus.

APPARATUS

Description	Quantity Required		Cat. No.
	Per Test	Unit	
Delivery Tubes, 90° with hook for TitraStir® Stir Plate	1.....	5/pkg.....	41578-00
Flask, Erlenmeyer, 125 mL.....	1.....	each.....	505-43
Flask, Erlenmeyer, 250 mL.....	1.....	each.....	505-46
Stir Bar, 28.6 x 7.9 mm.....	1.....	each.....	20953-52
TitraStir® Stir Plate, 115 Vac.....	1.....	each.....	19400-00
TitraStir® Stir Plate, 230 Vac.....	1.....	each.....	19400-10

5. To start titrant flowing and flush the delivery tube, hold the tip of the cartridge up. Advance the plunger release button to engage the piston with the cartridge (push the button in and toward the cartridge). Do not expel solution when pushing the piston toward the cartridge. Turn the delivery knob until air is expelled and several drops of solution flow from the tip. As you turn the knob a drive screw pushes a piston against the cartridge seal and forces liquid out through the delivery tube. Then use the counter reset knob to turn the digital counter back to zero and wipe the tip. The tip can be rinsed with deionized water rather than wiped, if desired.

GENERAL DESCRIPTION, continued

Figure 4 Using the TitraStir® Stir Plate





Figure 5 Using a Laboratory Stand



GENERAL DESCRIPTION, continued

Figure 6 Titrating the Sample



6. Use the smallest appropriate graduated cylinder or pipet to measure the sample volume from the given table. Transfer the sample into a 125-mL or 250-mL Erlenmeyer flask. Dilute to the appropriate total volume with deionized water if necessary.

Note: *Sample volume measurements and dilutions (if required) must be made accurately. However, final total volume of titrated solution is not critical.*

7. Add the necessary reagents to the sample and swirl to mix.
8. Immerse the delivery tube tip in the solution and swirl the flask while titrating. Titrate by turning the delivery knob. Keep turning the knob and swirling the sample until the end point is reached. Record the number of

GENERAL DESCRIPTION, continued

digits that appear in the digital counter window. See *Figure 6*.

Note: The number of digits required will usually range from 100 to 400. In nearly all of the procedures if the digits required is less than 100 or more than 400, an alternate sample volume or titrant cartridge should be used.

Note: Inaccurate results will occur if the delivery tube tip is held out of the solution rather than under the solution surface.

9. Calculate the concentration of your sample by using the following formula:

$$\text{Digits Required} \times \text{Digit Multiplier} = \text{Sample Concentration}$$

Where:

Digits Required = the number that appeared in the digital counter window in Step 8.

Digit Multiplier = the number from the table given in the procedure. It takes into account the sample dilution and titrant strength.

10. After completing testing for the day, press the plunger release button and manually retract the plunger into the body of the titrator. Remove the cartridge. Remove the delivery tube and reseal the cartridge with the polyethylene cap. See *Figure 7*.

Figure 7 Retracting the Plunger



11. Discard or clean the delivery tube immediately after use. To clean, force water, then air, into the tube opening with a syringe or wash bottle.

GENERAL DESCRIPTION, continued

1.3 Helpful Hints

1.3.1 To Reuse a Partially Emptied Cartridge

1. With the plunger fully retracted, attach cartridge to the titrator.
2. Press the plunger release; then manually push the plunger against the cartridge seal.
3. Attach a delivery tube. Hold the tip of the cartridge up. Eject air and a few drops of titrant, zero the counter, and wipe the tip.
4. Titrate as usual.

1.3.2 To Calculate Titrant Volume Used

Normalities of many Hach titration cartridge solutions have been designed so that the number of digits used in a titration corresponds to the sample concentration in mg/L. To determine the volume used in mL, divide the Digital Titrator reading by 800.

1.3.3 To Fill Your Own Titration Cartridges

Cartridges may be cleaned and refilled, or new empty cartridges, Cat. No. 14495-01, can be purchased from Hach Company. See *Figure 8*. When preparing to refill old cartridges, push the cartridge seal out of the cartridge with air pressure applied through the tip. Cap the tip, fill with solution and reinsert the cartridge seal using care to avoid wrinkling the Teflon sheath. Filling also can be accomplished at the tip with a syringe.

GENERAL DESCRIPTION, continued

Figure 8 Digital Titrator Cartridges



1.3.4 Verifying Technique

Whenever procedures are changed or new equipment is used, it is helpful to run a sample of known concentration. This technique will confirm the operator is following the procedure correctly and the new equipment is working properly. One objective important to Hach Company is making our tests self-verifying. This means Hach makes the tools available so the operator can check their own work for accurate results without relying on an outside lab or chemist.

For most of the tests in this manual, *Table 1* on page 21 lists each procedure, the suggested standard, the volume of standard needed, the titration cartridge used, and the number of expected digits when the test is performed correctly. The suggested standards are Voluette® or PourRite™ Ampules whenever possible because of their superior accuracy and stability.

To use titration standards follow these steps:

1. Select the procedure of interest and order the appropriate standard. Use the given catalog numbers.
2. Measure the volume of standard to be used as the sample in the procedure using a TenSette® Pipet or Class A pipet.
3. Perform the procedure as written, adding deionized water as necessary.

GENERAL DESCRIPTION, continued

- After titrating, the required number of digits should approximately equal the expected digits.

Call Hach Technical and Customer Service (1-800-227-4224) for additional help.

Table 1 Titration Standards

Procedure (Parameter)	Standard Description (Cat. No.)	Volume of Standard (mL)	Titration Cartridge (Cat. No.)	Expected Digits
Acid-Base: Acid	0.500 N H ₂ SO ₄ (2121-26)	1.0	1.600 N NaOH (14379-01)	250
		5.0	8.00 N NaOH (14381-01)	250
Base	0.500 N Na ₂ CO ₃ (14278-10)	1.0	1.600 N H ₂ SO ₄ (14389-01)	250
		5.0	8.00 N H ₂ SO ₄ (14391-01)	250
Acidity	0.500 N H ₂ SO ₄ (2121-26)	0.1	0.1600 N NaOH (14377-01)	250
		1.0	1.600 N NaOH (14379-01)	250
Alkalinity	0.500 N Na ₂ CO ₃ (14278-10)	0.1	0.1600 N H ₂ SO ₄ (14388-01)	250
		1.0	1.600 N H ₂ SO ₄ (14389-01)	250
Calcium*: mg/L CaCO ₃	10,000 mg/L CaCO ₃ (2187-10)	0.1	0.0800 M EDTA (14364-01)	100
		1.0	0.800 M EDTA (14399-01)	100
G.d.h.	10,000 mg/L CaCO ₃ (2187-10)	0.2	0.1428 M EDTA (14960-01)	112
		1.0	0.714 M EDTA (14959-01)	112
Carbon Dioxide	10,000 mg/L CO ₂ (14275-10)	0.2	0.3636 N NaOH (14378-01)	100
		2.0	3.636 N NaOH (14380-01)	

GENERAL DESCRIPTION, continued

Table 1 Titration Standards (Continued)

Procedure (Parameter)	Standard Description (Cat. No.)	Volume of Standard (mL)	Titration Cartridge (Cat. No.)	Expected Digits
Chloride	12,500 mg/L Cl (14250-10)	0.1	0.2256 N Hg(NO ₃) ₂ (14393-01)	125
		0.1	0.2256 N AgNO ₃ (14396-01)	125
		1.0	1.128 N AgNO ₃ (14397-01)	250
		1.0	2.256 N Hg(NO ₃) ₂ (921-01)	125
Chlorine	~50 mg/L Cl ₂ (14268-20) (see certificate)	2.0	0.02256 N Na ₂ S ₂ O ₃ (24091-01)	varies**
	~25 mg/L Cl ₂ (26300-20)	0.5	0.00564 N FEAS (22923-01)	varies***
Chromate	1000 mg/L Cr (2231 mg/L CrO ₄) (14664-42)	1.0	0.2068 N Na ₂ S ₂ O ₃ (22676-01)	223
Hardness: mg/L CaCO ₃	10,000 mg/L CaCO ₃ (2187-10)	0.1	0.0800 M EDTA (14364-01)	100
		0.1	0.0800 M CDTA (14402-01)	100
		1.0	0.800 M EDTA (14399-01)	100
		1.0	0.800 M CDTA (14403-01)	100
G.d.h.	10,000 mg/L CaCO ₃ (2187-10)	0.2	0.1428 M EDTA (14960-01)	112
		1.0	0.714 M EDTA (14959-01)	112
Iron	50 mg/L Fe (14254-10)	10.0	0.0716 M TitraVer (20817-01)	200
	1000 mg/L Fe (2271-42)	10.0	0.716 M TitraVer (20818-01)	100
Oxygen, Dissolved****	10 mg/L as DO (401-11)	100	0.2000 N Na ₂ S ₂ O ₃ (22675-01)	500
		200	2.00 N Na ₂ S ₂ O ₃ (14401-01)	100
Sulfite	5000 mg/L SO ₃ (22674-10)	1.0	0.3998 N KIO ₃ -KI (14961-01)	250

* One to two drops of Magnesium Standard Solution (10 g/L as CaCO₃) must be added to get a sharp end point. These added drops will not change the results.

GENERAL DESCRIPTION, continued

** The expected digits equal the volume of standard times the concentration on the certificate (e.g., 2 mL x 50 mg/L = 100 digits).

*** The expected digits equals the volume of standard times the concentration on the certificate times the constant, 4. (Example: 0.5 mL x 50 mg/L x 4 = 100 digits)

**** Add one Sulfamic Acid Powder Pillow to the volume of standard and follow Steps 10 to 12 in the Dissolved Oxygen Procedure. It is not necessary to add the first two reagents.

1.4 Adapting a Buret Titration to the Digital Titrator

Adapt any standard titration procedure using a buret to the Digital Titrator by using the following procedure.

1. Determine the approximate number of digits required. The Digital Titrator dispenses 1 mL per 800 digits on the counter. Using the following equation, determine the digits required for your buret method.

$$\text{Digits Required} = \frac{N_t \times mL_t \times 800}{N_c}$$

Where:

N_t = Normality of buret titrant

mL_t = milliliters of buret titrant required for an average titration

N_c = Normality of Digital Titrator cartridge

2. If the number of digits required is within the range of 70 to 350, you can use the procedure as written, substituting the Digital Titrator directly for the buret. Or, if the number of digits is outside of this range, make the following modifications:
 - a. If the number of digits required is more than 350, reduce the sample size to save titrant.
 - b. If the number of digits required is less than 70, increase the sample size to increase precision.
 - c. If the sample size is altered, adjust the amount of buffering or indicating reagents by the same proportion.
3. When using the Digital Titrator for your buret method, note the number of digits required for a sample titration.

GENERAL DESCRIPTION, continued

To convert the digits required to the equivalent number of milliliters if the buret method was used, calculate:

$$\text{Equivalent Buret Milliliters} = \text{Digits Required} \times \frac{N_c}{800 \times N_t}$$

If the sample size was changed, adjust the equivalent buret milliliters accordingly. If the sample size was increased, reduce the equivalent buret milliliters; if the sample size was reduced increase the equivalent buret milliliters. Multiply the equivalent buret milliliters by any normally used factors to calculate concentration in oz/gal, g/L, etc.

Example: Adapt a buret procedure, which normally requires about 20 mL of a 0.4 N titrant, to the Digital Titrator. Try an 8.0 N titration cartridge. The first equation above gives:

$$\text{Digits Required} = \frac{0.4 \times 20 \times 800}{8.0} = 800 \text{ digits}$$

Because this would use excessive titrant, reduce the sample size to one fourth its normal size to reduce the digits required to 200, well within the recommended range.


Upon completion of the titration using the smaller sample size, calculate the equivalent buret milliliters by the second equation above. If 205 were the digits required:

$$\text{Equivalent Buret Milliliters} = \frac{205 \times 8.0}{800 \times 0.4} = 5.13 \text{ mL}$$

Multiply the 5.13 mL by 4 to account for the reduction in sample size to give the true equivalent buret milliliters of 20.5 mL. If the buret method called for multiplying the number of milliliters of titrant by a factor to calculate the concentration of a sample component, then multiply 20.5 by that factor.

1.5 Using PermaChem® Powder Pillows

1. **Tap** the PermaChem on a hard surface to collect the powdered reagent in the bottom.




GENERAL DESCRIPTION, continued

2. **Tear** across on the dotted pillow line marked “TEAR” holding the pillow away from your face.



3. Using two hands, **Push** both sides toward each other until thumbs and forefingers form a diamond. Make sure to **Crease** the foil pack, so that it forms a spout.



4. **Pour** the pillow contents into the sample. The polyfilm lining is specially formulated to deliver all the powder necessary for accurate results (no tapping on the vessel edge is necessary).



1.6 Safety

Safety is the responsibility of each individual when performing analysis procedures, and the analyst must develop and maintain good safety habits. Because many of the procedures in this methods handbook use potentially hazardous chemicals and apparatus, it is important that the analyst practice good laboratory techniques to minimize accidents. The following paragraphs present several techniques applicable to water analysis in the laboratory and in the field. They are not all inclusive, of course, nor do

they apply only to the procedures provided in this handbook. They are general in nature but emphasize practices that are often key factors in personal injury incidents.

- Read labels carefully. Never remove the label from a reagent container. When preparing a reagent or standard solution, be sure to label the container clearly and date it.
- A Material Safety Data Sheet (MSDS) comes with each reagent. This sheet contains helpful information on first aid, spill and disposal procedures, and precautionary measures and should be read before using the product.
- Warning labels also appear on some of the apparatus used with the test procedures.
- Wear protective clothing when handling chemicals that cause irritation or burns. Eye protection in particular is important to guard against spattering and splashes from accidental spills when caustic materials are being used.
- Use tongs or finger cots when transferring apparatus that is hot.
- Use mechanical pipetters: Mouth pipetting could result in accidentally ingesting dangerous chemicals. Make a habit of using mechanical pipet fillers for all pipetting. This will avoid mistakes that could cause serious injury.
- Use special care with dangerous chemicals and apparatus.
- Follow the test procedure steps carefully and observe all precautionary measures. It is good practice to read the entire procedure carefully before beginning the procedure. Use safety equipment, such as pipet fillers, protective clothing, and ventilating hoods, appropriate for the test being conducted. Wipe up all spills promptly. Do not smoke or eat in an area where toxic or irritating chemicals are used. Use reagents and apparatus only as they were meant to be used and use them only as

GENERAL DESCRIPTION, continued

directed in the test procedure. Do not use damaged labware and malfunctioning equipment.

GENERAL DESCRIPTION, continued
