NL **LCK 334 Niet-ionische tensiden**

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GB **LCK 334 Nonionic surfactants**

! *Please check the "Edition Date" (see data table).*
■ *Safety advice and expiry date on package.*

Principle

Nonionic surfactants (ethoxylates with 3 to 20 ether bridges) react with cobalt thiocyanate, forming complexes. These are extracted in chloroform and photometrically evaluated.

Range of Application

Analysis of formulations and wash solutions; process analysis

Measuring range

Nonionic surfactants (sample = 1.0 mL)
Measuring range I.....**0.1 – 1.0 g/L**
Nonionic surfactants (sample = 0.2 mL)
Measuring range II.....**1 – 20 g/L**

Interferences

Cationic surfactants cause high-bias results to be obtained. High concentrations of anionic surfactants cause low-bias results to be obtained. Fatty acid esters and fatty acid alkanolamides are not detected. APGs (alkyl polyglycosides) are detected. Extremely slow phase separation and persistent streaks may make the analysis of highly polluted degreasing baths problematic or even impossible. In such cases the sample should be diluted appropriately (e.g. 1:100) and Cuvette Test LCK 333 should be used.

pH sample3 – 9
Due to the high salt content of acidic or alkaline degreasing baths, the pH must be adjusted to 7 – 8. Acidic degreasing baths with a high iron content should not be neutralized, as this would cause precipitation of iron hydroxide, which can absorb large amounts of surfactant. In such cases it is advisable to carry out a preliminary extraction with ethyl acetate, or to dilute the sample appropriately (1:100 - 1:500) and use the Cuvette Test LCK 333.

Temperature sample/reagents22°C
In case of not working at the right recommended temperature an incorrect result may be obtained.

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
Close cuvette and invert contents thoroughly.
Wait **2 min**.

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Kuвет **2 min** **krachtig schudden**.

Shake the cuvette **vigorously** for **2 min**.

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
NB: If streaks or droplets form in the lower part of the cuvette they can be eliminated by carefully tilting the cuvette back and forth and rotating it while holding it at an angle.

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
Then allow the cuvette to stand. Wait for the phase separation to finish (**a few minutes**).

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Thoroughly clean the outside of the cuvette and evaluate.

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Note:

Depending on the surfactant types and their solubility, the phase separation may proceed very slowly or not at all. In such a case, the addition of 0.3 mL ethylene glycol (dealer in laboratory products) to the cuvette immediately before adding the sample may accelerate the phase separation. If this method is used, it is advisable to carry out a specific calibration with the surfactants.

	Analysenküvette ① Cuve d'analyse Cuvetta d'analisi Analyse-kuvet Sample cuvette
 Barcode ¹⁾	✓

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